



Test Report Summary of SR SMITH Pool Lift Model No. AXS1000L

To SR Smith LLC

Attention: Mr. Bill Svendsen

Address: 1017 SW Berg Parkway
P.O. Box 400
Canby, OR 97013

Email: BSvendsen@srsmith.com

Reference No.: MEIC-7805006-SRS

Date: August 20, 2012

**Subject: Summary Report of Swimming Pool Lift Tests per ADA 1009.2
Model No. AXS 100L**

This summary contains the limited results of the detailed report that was prepared as an outcome of the testing performed on a sample of swimming pool lift **SR Smith Model No. AXS1000L, AXS Aquatic Lift**, which was tested according to the American's with Disabilities Act (ADA) standard 1009.2. All the tests were performed under simulated pool environment and on site at SR Smith, LLC facility, located in Canby, Oregon.

Scope of the Testing (Work)

The testing was conducted under simulated pool environment to ascertain that the pool lift (tested in as-mounted condition by SR Smith) complied with ADA Standard 1009.2.

Mounting Apparatus

A platform made of plywood simulating a typical pool deck was set up by SR Smith, and was used to mount the pool lift. The pool lift to be tested was mounted into a simulated concrete instillation at a typical installation distance (18" on center) to the simulated pool wall.

Tests Conducted and Measurements Made

A test protocol to include all the measurements and tests as detailed in the standard was prepared and tests were performed by strictly following the Standard protocol and by using calibrated instruments. All the measurements (where applicable) were performed by using the deck to mimic the pool, and the lifts mounted on it.

Section 1009.2.1 – Scope of the test: Pool Lift Location (not conducted)

Section 1009.2.2 - Scope of Test: Seat Location.

Section 1009.2.3 - Scope of Test: Clear Deck Space.

Section 1009.2.4 - Scope of Test: Seat Height.

Section 1009.2.5 - Scope of Test: Seat Width.

Section 1009.2.6 - Scope of test: Footrests and Armrests.

Section 1009.2.7 - Scope of Test: Operation.

Section 1009.2.8 - Scope of Test: Submerged Depth

Section 1009.2.9 - Scope of Test: Lifting Capacity

MEI-Charlton, Inc.

Website: www.meic.com

7220 N Lombard St, Portland, OR 97203-3208

Tel.: 503-228-9663; Fax: 503-228-4065

TO: SR Smith LLC
 SUBJECT: Pool Lift Tests per ADA 1009.2
 REF. NO.: MEIC-7805006-SRS

Table 1: Test Results Data of Swimming Pool Lift Tests per ADA 1009.2, for, SR Smith Model No. AXS 1000L,

ADA Req.	Test Details	Observation	Req.	Pass	Fail
1009.2.1	Pool Depth at Installation Location	***	≤48 inches		
1009.2.2	Edge of Pool to Centerline of Seat Distance when Raised	18 inches	≥16 inches	X	
1009.2.2	Deck Slope-Edge of Pool to Centerline	***	≤1:48		
1009.2.3	Clear Deck Space Width	-	≥36 inches	X	
1009.2.3	Clear Deck Space Length from 12 in Behind Rear of Seat	-	≥48 inches	X	
1009.2.3	Deck Slope-Clear Deck Space	***	≤1:48		
1009.2.4	Seat Height from Deck to Top of Seat Surface when Raised	21 ¾ inches maximum	16-19 inches	X	
1009.2.5	Seat Width	18 ½ inches	≥16 inches	X	
1009.2.6	Footrests Provided (y/n)	Yes	Yes	X	
1009.2.6	Footrests Move with Seat (y/n)	Yes	Yes	X	
1009.2.6	Armrests Provided (y/n)	-	NR		
1009.2.6	Armrest on Side Opposite Water Removable (y/n)	-	Yes*		
1009.2.6	Armrest on Side Opposite Water Folds Clear of the Seat (y/n)	-	Yes*		
1009.2.7	Pool Lift Operable from Deck Level (y/n)	Yes	Yes	X	
1009.2.7	Pool Lift Operable from Water Level (y/n)	Yes	Yes	X	
1009.2.7	Controls and Operating Mechanism Unobstructed When Lift is in Use (y/n)	Yes	Yes	X	
309.4	Operable Parts Operable with One Hand (y/n)	Yes	Yes	X	
309.4	Tight Grasping Required (y/n)	No	No	X	
309.4	Pinching Required (y/n)	No	No	X	
309.4	Twisting of the Wrist Required (y/n)	No	No	X	
309.4	Force Required to Activate Button Less than 5 lbs (y/n)	0.75lb 0.75lb 1.20lb 1.20lb	≤5lb	X	
1009.2.8	Depth of Seat Submersion	19 ¼ inches	≥18 in	X	
1009.2.8	Depth of Seat From Deck	25 ¼ inches	≥18 in	X	
1009.2.8	Edge of Deck to Centerline of Seat Distance when Lowered	***	NA		
1009.2.9	Lift Sustains Static Loading of 450 lbs or 1.5 x Rating (y/n)	450lb	Yes	X	

TO: SR Smith LLC
SUBJECT: Pool Lift Tests per ADA 1009.2
REF. NO.: MEIC-7805006-SRS

Page 3 of 3

Notes: Rated at 300lbs.

No. *If armrests are provided they must be either removable or able to fold clear on side opposite water.

**Armrest is removable by removing screws. ; NA – Not Applicable; NR – Not Required

***Item being a pool specific requirement was not tested.

+ Req = Requirement

Results: Swimming Pool Lift, Instrument #7, SR Smith Model No. AXS1000, SR Smith AXS Aquatic Lift, in the simulated pool environment was found to pass all tests that were performed.

Comment: SR Smith Sample Model No AXS1000 passed the tests as identified in the above table.

This summary report is based on the full report which contains the details of the tests conducted on the identified model and as supplied by SR SMITH, LLC presents the outcome of the testing performed by MEIC. This summary is not a substitute for the full report.

A copy of the full report is available to all SR SMITH Prospective and Current customers from MEIC at a nominal cost. You may order your copy at www.meic.com directly.

Disclaimer:

© Copyright 2012, MEI-Charlton, Inc (MEIC). This report constitutes a confidential communication between the client and MEIC, and the report, as a whole or in parts is not allowed be published in any form or distributed to any one without MEI-Charlton, Inc's written permission. This report is meant for information purpose only, and to inform the identified client the outcome of the study/testing commissioned by the client explicitly identified at the beginning of the report. Any questions, issues related to the scope of work carried out by MEIC and reported in this report must be conveyed in writing to MEIC within 30 days of the date of issue of this report (sent to client by electronic or other means). Absence of such issues, questions will be an explicit acceptance of the report in as is form. All additional consultation time, experiments will be in addition to the original quote and invoiced amount and will be charged at the prevailing rates. Client and/or their authorized representatives may choose to use the results of this report at their own risk and have agreed to hold MEIC, its officers, employees, principals, stock holders (past, present and future) and their successors free of any responsibility whatsoever. Client further agrees to reimburse any damages and the costs fees to defend that MEIC may incur if a law suit or damages are awarded against MEIC, hereunder any prevailing laws.

MEI-Charlton, Inc.

Website: www.meic.com

7220 N Lombard St, Portland, OR 97203-3208

Tel.: 503-228-9663; Fax: 503-228-4065